




Marker Compatibility for Stereotactic Breast Biopsy



Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape	Metal	Bio-absorbable Carrier	Deployment Method	# per box	
Brevera	TUMARK-BREV-P-VISION	18G	12.3cm	VISION	Nitinol	No	Introducer	10	
	TUMARK-BREV-P-X			X	Nitinol	No	Introducer	10	
	TUMARK-BREV-P-Q		12.5cm	Q	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-VISION		12.7cm	VISION	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-X			X	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-Q		12.9cm	Q	Nitinol	No	Introducer	10	
Eviva and Brevera	Smark-Eviva-10	12G	10cm	Mini Cork	Ti	Yes	Introducer	10	
	Smark-Eviva-2s-10			Top Hat	Ti	Yes	Introducer	10	
	Smark-E13-ss1		13cm	Buckle	SS	Yes	Introducer	10	
	Smark-E13-ss2			Infinity	SS	Yes	Introducer	10	
	Smark-E13-ss3			Stoplight	SS	Yes	Introducer	10	
	Smark-Eviva-13			Mini Cork	Ti	Yes	Introducer	10	
	Smark-Eviva-2s-13		Top Hat	Ti	Yes	Introducer	10		
	TriMark-Eviva-10		15G	10cm	Cork	Ti	No	Introducer	10
	TriMark-Eviva-2S-10	Hourglass			Ti	No	Introducer	10	
	TriMark-Eviva-13	13cm		Cork	Ti	No	Introducer	10	
	TriMark-Eviva-2S-13			Hourglass	Ti	No	Introducer	10	
	Eviva	TUMARK-E13-S-VISION	18G	12.9cm	VISION	Nitinol	No	Introducer	10
		TUMARK-E13-S-X			X	Nitinol	No	Introducer	10
		TUMARK-E13-S-Q		13.2cm	Q	Nitinol	No	Introducer	10
TUMARK-E13-P-VISION		13.3cm		VISION	Nitinol	No	Introducer	10	
TUMARK-E13-P-X				X	Nitinol	No	Introducer	10	
TUMARK-E13-P-Q		13.5cm		Q	Nitinol	No	Introducer	10	

Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape	Metal	Bio-absorbable Carrier	Deployment Method	# per box
ATEC	TriMarkTD-2S-13-09	9G	13cm	Hourglass	Ti	No	Needle Hub	10
	Smark-ATEC-13-09			Mini Cork	Ti	Yes	Needle Hub	10
	Smark-ATEC-2s-13-09			Top Hat	Ti	Yes	Needle Hub	10
	Smark-A3609-ss1		36cm	Buckle	SS	Yes	Handpiece	10
	TriMark TD 1309			Cork	Ti	No	Needle Hub	10
	Smark-Atec-2s-36-09			Top Hat	Ti	Yes	Handpiece	10
	Smark-ATEC-36-09			Mini Cork	Ti	Yes	Handpiece	10
	TriMark TD 3609			Cork	Ti	No	Handpiece	10
	TriMarkTD-2S-36-09			Hourglass	Ti	No	Handpiece	10
	351252			38.2cm	U	Nitinol	No	Handpiece
	TriMark TD 1312	12G	13cm	Infinity	SS	Yes	Needle Hub	10
	TriMarkTD-2S-13-12			Cork	Ti	No	Needle Hub	10
	Smark-A1312-ss2			Hourglass	Ti	No	Needle Hub	10
	Smark-ATEC-13-12			Mini Cork	Ti	Yes	Needle Hub	10
	Smark-ATEC-2s-13-12			Top Hat	Ti	Yes	Needle Hub	10
	Smark-ATEC-36-12		36cm	Mini Cork	Ti	Yes	Handpiece	10
	TriMark TD 3612			Cork	Ti	No	Handpiece	10
	TriMarkTD-2S-36-12			Hourglass	Ti	No	Handpiece	10
	351250			38.2cm	U	Nitinol	No	Handpiece

		
<p>Tumark® Marker Intelligently designed to provide long term ultrasound visibility, non-bio-absorbable, bio-compatible permanent marker with excellent visibility under mammography at deployment.¹</p>	<p>SecurMark® Marker The bio-absorbable suture-like net that carries a bio-compatible permanent marker, is designed to minimize movement within the biopsy cavity and enhance ultrasound visualization for three to four weeks post-biopsy.²</p>	<p>TriMark® Marker Reliable and simple bio-compatible Titanium metal markers, optimal for thin-breasted patients and superficial lesions.</p>

1. In 85% of marker placements, physicians stated that the ultrasound visibility was good to excellent upon deployment. In 56% of marker placements, physicians stated that ultrasound visibility was excellent upon deployment. Based on Tumark® Marker Data Collection Study DHM-06169, 3 clinicians at 3 hospitals for 90 marker placements. 2. Based on Clinical Evaluation Report DHM-08585. Visibility is dependent on surrounding tissue, experience may vary.

www.hologic.com | BSHInsideSales@hologic.com | +1.800.442.9892

Hologic Headquarters
United States / Latin America
250 Campus Drive
Marlborough, MA 01752
USA Tel: +1.508.263.2900
Sales: +1.781.999.7453
Fax: +1.781.280.0668
Email: info@hologic.com

Hologic Europe
Everest (Cross Point)
Leuvensesteenweg
250A 1800
Vilvoorde, Belgium
Tel: +32.2.711.4680
Fax: +32.2.725.2087

Hologic Asia Pacific Limited
Unit 01-03A, 13/F
No. 909 Cheung Sha Wan
Road Cheung Sha Wan
Kowloon, Hong Kong
Tel: +852 3748 7700
Fax: +852 3526 0723

Hologic Australia
Suite 402, Level 4
2 Lyon Park Road
Macquarie Park NSW 2113
Australia
Tel: +61.2.9888.8000
Fax: +61.2.9870.7555

MISC-01695 Rev.010 (4/20) ©2020 Hologic Inc. All rights reserved. Hologic, ATEC, Eviva, Brevera, SecurMark, TriMark, Tumark and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.